

Technical datasheet

Object: "Saliva Ejector (for dental use)"

Use: these products are used exclusively by Dentists or qualified personnel in Dental Surgeries during the various operations needed to extract saliva from the oral cavity of the patient

Material: they are small PVC tubes with a core of copper and a perforated cap inserted at one extremity of the tube (the other end of the tube is connected to an aspirator system).

Size: length 13 cm and 14,5 cm, diameter 6,5 mm.

Cap length: 15 mm

Characteristic: the product is totally latex free, glue free, with no toxic substances or contamination. The characteristic of this product is that the copper core of the cannula permits the Dentist or qualified personnel to obtain the form required for the patient under exam by simply manipulating with their hands. With smooth tip and shaft, tip is securely bonded to shaft.



Fig.1 saliva ejector before manipulating



Fig.2 saliva after manipulation

Classification: Class: II A medical device according to rule 5 (last paragraph) of the legislative decree 46/97, in acknowledgement of the European Directive 93/42/CEE and successive additional modifications – ref. European Directive 2007/47/CE.

Packaging: 100pcs. in sealed polyethylene packets.

Method of Conservation: Away from direct heat

Chiasso, 15.02.2014

AKZENTA INC.
Branch Lugano
Dr. Joachim Missfeldt